

# E12.5 & E14.5 mouse DRG dissociation for FACS, 053114

Masato Hoshi, MD, PhD,  
Sanjay Jain Lab, Renal Division, Department of Medicine,  
Washington University School of Medicine

## Reagents:

- Hibernate-E; Life Technologies, #A12476-01
- DMEM/F12; Life Technologies, #11330
- 0.05% Trypsin/0.02% EDTA; Sigma, #59417C or #T3924
- DNaseI; Sigma, #D5025
- RNase inhibitor; RNaseOUT (Life Technologies, #10777, 40units/ul) or equal
- 7-AAD; Life Technologies, #A1310 (powder, dissolve in DMSO at 1mg/ml)

**Reference: modified from** (Sasaki et al, J. neuro 2009).

## Procedures:

- 1) Euthanize a pregnant mouse following your animal protocol approval and extract embryos/fetuses from the uterus.
- 2) Save the embryos/fetuses in a 24-well plate with 1ml of Hibernate-E medium (Life Technologies) in each well of the plate on ice.
- 3) Take one embryo in 10cm plastic plate with dissection medium (DMEM/F12 (Life Technologies) or Hibernate-E (30ml)) and dissect the DRGs from embryos. (T11 to S1, total 20 DRGs from each embryo; or as needed for your study)  
(\*when you use DMEM/F12, the medium should be put in a CO<sub>2</sub> incubator (cell culture, 5%CO<sub>2</sub> incubator at 37C is fine) for more than 2 hours to keep the pH balanced (gauge by the color, prefer orange).
- 4) Collect DRGs in a 1.5ml tube on ice with dissection medium in a pipette tip. Keep the tube on ice until all dissections are done.
- 5) Centrifuge at 500 xg for 5min at 4C and discard as much supernatant as possible by using P-200 pipette.
- 6) Add 500ul of 0.05% Trypsin/ 0.02% EDTA (Sigma) with 200ug/ml of DNaseI (Sigma) to the tube and briefly vortex to mix.
- 7) Incubate it at 37C for 10min.
- 8) Triturate with P-1000, then P-200 pipette till you notice complete dissociation of the specimen. Add 500ul of DMEM/F12 + 10% FCS (filtered) and mix by pipetting.
- 9) Centrifuge at 500 xg for 5min at 4C and discard supernatant by using a pipette. Leave around 30ul to avoid disturbing the cell pellet.
- 10) Add 500ul of PBS + 5% FCS (filtered) with 0.1units/ul of RNase inhibitor (Life Technologies, 40units/ul, 25ul to 10ml PBS/FCS solution) and mix by pipetting.
- 11) Centrifuge at 500 xg for 5min at 4C and discard supernatant by using a pipette.
- 12) Repeat this washing step twice more with PBS/FCS/RNase inhibitor solution. Total three times. In the last wash, the cells get dispersed very quickly with only few rounds of pipetting.
- 13) After discarding supernatant, add 300ul (or appropriate volume for your FACS) of PBS + 5% FCS (filtered) with 0.1units/ul of RNase inhibitor and mix well by pipetting.

- 14) Filter the dissociated cells with 40um pore size cell strainer (BD, #352340) and collect the cells into a FACS tube (BD, #352063).
- 15) Add 7-AAD at 1ug/ml concentration, if needed, and do sorting. 7-AAD is added 10min before sorting. We have had good success using MoFlo, Sony or Ariall, at 100 um nozzle size, and 30psi.
- 16) Sorted cells are collected in 750ul of TRIzol-LS (Life Technologies, #10296-010) for RNA extraction and saved in -80C until use.